

PRESTO™ A80t

Heating a 20 liters reactor from 0 °C to +20 °C

Objective

This case study tests the heating power of PRESTO™ A80t with a 20 liters glass reactor. The PRESTO™ A80t is connected to the reactor via two 1 m metal tubings. The PRESTO™ A80t is programmed to heat up from 0 °C to +20 °C.

Environment

Room temperature +20 °C
Humidity 45%
Voltage 208 V / 60 Hz

Test Conditions


JULABO unit	PRESTO™ A80t
Cooling power	+20 °C 1.2 kW
	0 °C 1.2 kW
	-20 °C 1.1 kW
Heating capacity	3.4 kW
Band limit	with
Flow pressure	0.5 bar
Bath fluid	Thermal HL80
Reactor	20 liters glass reactor (Asahi) filled with 19 l Thermal HL80
Jacket volume	7 l
Control	External (ICC)



Test Results

The PRESTO™ A80t heating process from 0 °C to +20°C in 27 min without overshoot.

Measured with *EasyTEMP Professional*



- Setpoint
- Temperature in reactor's interior
- Temperature in reactor's jacket

Tip

You can also use the robust Pt100 with PTFE coating.



Tip

Use our tube adapters and your tubing will no longer kink.

