

Julabo Case Study

JULABO PRESTO® A40

Heating a 5 liters reactor
from +100 °C to +200 °C



Objective


This case study tests the heating power of JULABO PRESTO® A40 with a 5 liters glass reactor. The A40 is connected to the reactor via two 2.0 m metal tubings. The A40 is programmed to heat up from +100 °C to +200 °C.

Test Conditions

| | |
|------------------|--|
| JULABO unit | JULABO PRESTO® A40 |
| Cooling power | +20 °C 1.2 kW |
| | 0 °C 0.9 kW |
| | -20 °C 0.6 kW |
| Heating capacity | 2.7 kW |
| Band limit | No |
| Flow pressure | 0.40 bar |
| Bath fluid | JULABO Thermal HL40 |
| Reactor | 5 liters glass reactor (Rettberg) filled with 5 liter JULABO Thermal HL40 |
| Control | External (ICC) |

Environment

| | |
|------------------|---------------|
| Room temperature | +20 °C |
| Humidity | 45 % |
| Voltage | 230 V / 50 Hz |



Test Results

See chart on back page: The A40 heating process from +100 °C to +200 °C in 45 min without overshoot.


Tip

You can also use the robust Pt100 with PTFE coating.

More tips on back page >>



JULABO GmbH
Eisenbahnstraße 45
77960 Seelbach / Germany
Tel. +49 (0) 7823 51-0



- Setpoint
- Temperature in reactor's interior
- Temperature in reactor's jacket

Tip

Make use of the option to regulate the pump pressure. You can define the desired pressure in the PRESTO® settings.



Tip

The Ethernet interface permits full access to all operational functions of the PRESTO®.



JULABO GmbH
 Eisenbahnstraße 45
 77960 Seelbach / Germany
 Tel. +49 (0) 7823 51-0